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A handwritten signature in blue ink, appearing to be 'JGA'.

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[www.gosdravnadzor.ru](http://www.gosdravnadzor.ru)

**CLINICAL EVALUATION REPORT**

<b>INTERNAL PROJECT NUMBER:</b> BTL-049		<b>PRODUCT BUSINESS NAME:</b> BTL-6000 TR-Therapy	
<b>Version:</b> 103		<b>Signature:</b>	<b>Date:</b>
<b>Made by:</b> <b>Title:</b>	Zuzana Heinzová		4.11.2014
<b>Reviewed by:</b> <b>Title:</b>	Lukáš Množil		4.11.2014
<b>Approved by:</b> <b>Title:</b>	Ondřej Prouza		4.11.2014

**RECORD OF CHANGES**

Revision	Date	Last version			Change Description
		Version	Valid from	Valid to	
1	October 24, 2014	100	October 1, 2014	October 24, 2014	The 2 product variants were specified. The variants differ only by the display dimension. This difference doesn't affect the clinical performance and safety of the product
2	October 30, 2014	101	October 24, 2014	October 30, 2014	The revision of articles and equivalent devices was made.
3	November 4, 2014	102	October 30, 2014	November 4, 2014	Update due to traceability purposes

Director BTL Industries Limited:

\_\_\_\_\_/Ondřej Sojka

**BTL Industries Limited**  
 161 Cleveland Way  
 Stevenage  
 SG1 6BU Hertfordshire  
 United Kingdom

220

Clinical Evaluation Report  
BTL-6000 TR-Therapy

Content

1	Introduction.....	4
2	General Description of BTL-6000 TR-Therapy and its intended use.....	4
2.1	Device Model.....	4
2.2	Device Classification.....	4
2.3	Device Description.....	5
2.3.1	Front Panel of the BTL-6000 TR-Therapy.....	5
2.3.2	Rear Panel of the Device.....	5
2.3.3	Applicators and electrodes.....	7
2.3.3.1	Neutral electrode.....	7
2.3.3.2	Applicator with the active electrode.....	7
2.3.3.3	Active electrode.....	7
2.3.4	Active electrode sizes.....	7
2.4	Intended use and indications.....	8
2.5	Intended population.....	8
2.6	Contraindications.....	8
2.7	Possible side effects.....	8
3	Clinical Evaluation Process.....	9
3.1	Clinical Evaluation Justification.....	9
3.2	Clinical Evaluation Objectives.....	9
4	Equivalent devices.....	10
4.1	Comparison with equivalent device.....	10
4.2	Clinical Experience with equivalent devices.....	10
5	Clinical safety.....	12
5.1	Equivalent device.....	12
5.2	Post Market Surveillance Data.....	12
5.3	Risk Analysis Review.....	12
5.3.1	Residual risks.....	12
5.3.2	Clinical risks.....	13
5.4	Conclusion – clinical safety.....	13
6	Clinical performance.....	14
6.1	Equivalent device.....	14
6.2	List of components.....	14
6.3	List of internal and external tests.....	14
6.4	Manufacturing process.....	14
6.5	Conclusion – clinical performance.....	15
7	User's manual.....	16
7.1	User's manual - Device intended use and contraindication.....	16
7.2	User's manual - Warnings.....	16

Handwritten initials

Clinical Evaluation Report  
BTL-6000 TR-Therapy

7.3	Conclusion – User’s manual .....	16
8	Conclusion.....	17
8.1	Intended Use Benefit evaluation .....	17
8.2	Conclusion - Clinical risk / benefit .....	17
8.3	Overall evaluation .....	17
9	Documentation – cross reference .....	18
9.1	Unpublished literature .....	18
10	Statement.....	18

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Clinical Evaluation Report  
BTL-6000 TR-Therapy

**1 Introduction**

The *Clinical Evaluation Report* is compiled for the purpose of providing the evidence of effective and safe clinical use of BTL-6000 TR-Therapy.

**2 General Description of BTL-6000 TR-Therapy and its intended use**

BTL-6000 TR-Therapy is a professional physiotherapy device. It is equipped with a colour touch screen, which considerably simplifies its operation. The on-screen information will guide the user step-by-step through the entire therapy process. The therapeutic parameters are easily set using the touch screen buttons and knobs/keys on the device.

Therapy is simply started by manually set using the touch screen buttons. Throughout the therapy the device keeps you informed on the screen about the applied therapy type, remaining therapy time and main therapy parameters.

BTL-6000 TR-Therapy is device for tissue heating by means of a 500 kHz contact radiofrequency. The therapist works with device generating high-frequency current. One of the electrodes, plate shaped, is neutral (passive), therapist places the electrode under one part of the patient's body. Second electrode is the active electrode and performs the therapy.

The device is available in two variants BTL-6000 TR-Therapy Elite and BTL-6000 TR-Therapy Pro. The difference between the two variants is in the size of the display. BTL-6000 TR-Therapy Elite has 8.4" touch screen and BTL-6000 TR-Therapy Pro has 5.7" touch screen.

**2.1 Device Model**

<b>Product</b>	BTL-6000 TR-Therapy
<b>Models</b>	BTL-6000 TR-Therapy Elite BTL-6000 TR-Therapy Pro

**2.2 Device Classification**

The BTL-6000 TR-Therapy is an active medical device. It is intended to control the performance of an active therapeutic device in Class IIa.

As such, it falls into Class IIa, pursuant to the Rule 9 of Annex IX of the MDD in effective version.

Pursuant to Article 11, Section 2, second paragraph of the MDD, the Device follows Annex II route of conformity, i.e. full quality assurance (point 4 of Annex II is not applicable).

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
113

Clinical Evaluation Report  
BTL-6000 TR-Therapy

2.3 Device Description

2.3.1 Front Panel of the BTL-6000 TR-Therapy


1. The touchscreen (8.4 " for BTL-6000 TR-Therapy Elite, 5.7 " for BTL-6000 TR-Therapy Pro)
2. The **select** knob: To select individual parameters
3. The **enter** key
4. The **esc** key
5. The **start/stop** key: To start and stop therapy
6. The **on/off** switch: Back-lit in blue when the control unit is switched on
7. The USB port: Located in the space used for gripping and lifting the device. The USB port is used ONLY for technician service purposes such as the uploading of firmware. It is not designed for therapy use!

2.3.2 Rear Panel of the Device

114

Clinical Evaluation Report  
BTL-6000 TR-Therapy

- 8. Neutral electrode output connector (black)
- 9. Capacitive applicator output connector (yellow)
- 10. Resistive applicator output connector (blue)
- 11. Mains power switch
- 12. Mains cable connector
- 13. Fuse
- 14. Fuse

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
Handwritten signature or initials.

Clinical Evaluation Report  
 BTL-6000 TR-Therapy

### 2.3.3 Applicators and electrodes

There are two possible types of applicators. Each applicator has on its end a sliding type of connector. Device uses two types of electrodes: active electrodes and neutral electrodes. Active electrodes can be either capacitive or resistive, each of them in three different sizes. The active electrodes are used for the treatment itself. The neutral electrode is used as a reference electrode.


#### 2.3.3.1 Neutral electrode

#### 2.3.3.2 Applicator with the active electrode

#### 2.3.3.3 Active electrode

### 2.3.4 Active electrode sizes

Electrodes	Capacitive	Resistive
Size	. 45 mm . 65 mm . 80 mm	. 45 mm . 65 mm . 80 mm

## Clinical Evaluation Report BTL-6000 TR-Therapy

### 2.4 Intended use and indications

The device is intended for treatment of musculoskeletal system acute, chronic and degenerative disorders and injuries such as: sports pathologies of bones, joints and muscles; tendinous, muscular and articular pathologies; neck pain; pain in patients with isolated long head biceps tendinopathy; post operative pain and oedema treatment of femur fractures.

### 2.5 Intended population

The device is intended for the adult patients.

### 2.6 Contraindications

Contraindications for contact radiofrequency current: Patients with a pacemaker or other electronic implants must be approved beforehand by a specialist in a particular medical field. Patients must not touch the electrical conductors or grounded object, formation of undesirable pathways could occur for the transmission of high-frequency currents. People with low blood pressure or patients with treatment of the cervical vertebrae, spine or joints have to take a longer lasting therapeutic procedure, the therapy can lead to further reduction of blood pressure. Use of the device is contraindicated for pregnant women. When using resistance mode, it is necessary to remove the dentures if treatment is carried out in the relevant area. Patients or the users should remove rings, necklaces, piercings and other metal objects. Do not apply on damaged skin, for example in the case of open wounds or fresh burns.

### 2.7 Possible side effects

Possible side effects for contact radiofrequency current: Unpleasant sensation can temporarily occur in the treated area, moderate skin irritation.

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228

## Clinical Evaluation Report BTL-6000 TR-Therapy

### 3 Clinical Evaluation Process

In accordance with the definition of clinical data set forth in Medical Devices Directive 93/42/EEC, as amended by 2007/47/EC, we have chosen the way of collecting clinical information about safety and performance following from a device's intended use: Through published and/or unpublished written documents concerning clinical experience with equivalent devices. The literature search (*Comparative literature study for BTL-6000 TR-Therapy*) was made to find out the devices placed on the medical device market for a long time which have the equivalent parameters as the BTL-6000 TR-Therapy and the data supported their intended use. The chosen devices and the articles supported particular indication from their intended use are further evaluated within this clinical evaluation report. Four equivalent devices were found. For each of them we found the articles which support particular intended use of the devices. The particular indication from the articles was chosen and based on them we set the intended use for the BTL-6000 TR-Therapy.

#### 3.1 Clinical Evaluation Justification

We reviewed the actual published literature within the literature search for the up to date information about the devices offers a non-invasive therapeutic method based on a conversion of a radiofrequency current to heat energy. Within the literature search we found the equivalent devices which we used as a comparable devices to based the clinical evaluation on. Within the literature search we didn't find any new information and any new clinical risk then is stated for the equivalent devices. So, we could base the clinical evaluation of BTL-6000 TR-Therapy on the long term clinical experience with devices which have the same technology principles and same intended use. The long term manufacturing and clinical experience of Company, manufacturer the tecar therapy equivalent devices, from the market allow us to use the data to demonstrate the clinical performance and safety of the evaluate device. Performing of clinical study was not required, because of sufficient data obtained within the literature search which we can use to based the intended use on.

#### 3.2 Clinical Evaluation Objectives

In accordance with the definition of clinical data set forth in Medical Devices Directive 93/42/EEC, as amended by 2007/47/EC, we have chosen the way of collecting clinical information about safety and performance following from a device's intended use: Through published and/or unpublished written documents concerning clinical experience with equivalent devices manufactured by Indiba S.A. The clinical evaluation is based on clinical, technical and biological equivalency proved for the equivalent devices and BTL-6000 TR-Therapy. In accordance with the definition of clinical data set forth in Medical Devices Directive 93/42/EEC, as amended by 2007/47/EC, we have chosen two ways of collecting clinical information about safety and performance following from a device's intended use: Through published and/or unpublished written documents concerning clinical experience either with an equivalent competitive devices for which clinical, technical and biological equivalence to the evaluated device can be demonstrated. The critical evaluation of currently available relevant scientific literature was made to find the relevant clinical data.

The entire process of how to select equivalent devices and a literature search are described in the *Comparative Literature Study for the BTL-6000 TR-Therapy*.

For evaluation of comparable clinical data such equivalent devices have been used, which show the highest number of identical or possibly similar parameters to the Device. The analogy/equivalence of the BTL-6000 TR-Therapy device to currently used equivalent devices is based on clinical, biological, technical equivalence.

228

## Clinical Evaluation Report BTL-6000 TR-Therapy

### 4 Equivalent devices

BTL-6000 TR-Therapy has many competitive devices being on the market for many years as well. The manufacturer of all devices found out within the literature search is the Spain origin company Indiba S.A. The Italian distributor for this manufacturer is named Human Tecar. This is why some of the products can be found as Human Tecar company products and some can be found as Indiba company product. In the studies are mentioned following 4 products: Human Tecar HCR 900, 901, TECAR 901 and MD-308.

Since 1983, **INDIBA**® is a company that specializes in research, development, manufacture and global distribution of medical and aesthetic radiofrequency devices based on a unique exclusive and patented technology. This electro regenerative technology provides beneficial biological effects in tissue regeneration, pain relief and cicatrisation. The latest research of the company is focused on ion exchange and cellular electrical potentials, giving rise to the Proionic® System. Clinical references using our aesthetic and therapeutic radiofrequency, demonstrating the efficiency of the system around **INDIBA**® treatment protocols, are supported by more than 150 scientific publications in the fields of plastic surgery, aesthetic medicine, sports medicine, traumatology and physiotherapy.

Based in Barcelona and with over 30 years' experience, we are delighted that **INDIBA**® has made its treatments available worldwide with over 15,000 units of devices installed. **INDIBA**® has been present in Italy for 25 years, in Japan for 20 years and is successfully pursuing a targeted expansion across Europe, Asia and Latin America. One of our core values is to provide our partners and clients with a unique technology that consistently provides effective and visible results. In order to underpin our success we always provide our professional customers with a first class education and customized training program in order to achieve total satisfaction from our patients and clients.

The tecar therapy device stimulates the natural mechanisms of healing and tissue recovery non-invasively. It uses a specific form of high frequency current ( $490 \pm 5\%$  kHz) applied with a system combining the capacitive and resistive modes. It speeds up regeneration, scar formation and has analgesic, anti-inflammatory and fibrolytic properties. The radiofrequency current is transferred to a human body through the electrodes connected to the device. In the human body the radiofrequency current is converted to the heat energy, which influences the tissue.

#### 4.1 Comparison with equivalent device

The analogy/equivalence of the BTL-6000 TR-Therapy device to currently used equivalent devices is based on:

##### Clinical equivalence

The BTL-6000 TR-Therapy is used for the equivalent clinical condition and purpose as HCR 900, HCR 901, TECAR 901, MD-308.

##### Biological equivalence

All devices used the same type of applicators.

##### Technical equivalence

The BTL-6000 TR-Therapy works with the same technical principles and comprises the equivalent technical components as the HCR 900, HCR 901, TECAR 901, MD-308.

The comparison is described in the *Comparative Literature Study for BTL-6000 TR-Therapy* in the chapter 3.5.2

#### 4.2 Clinical Experience with equivalent devices

The main goal of the clinical literature search is fulfilled. The equivalent devices placed on the market for a long time were found. The clinical data supported their intended use are shown in the obtained whole articles.

The equivalency is proved based on the clinical, biological and technical comparison. We can declare regarding this comparison of technical, biological and clinical parameters described above that BTL-6000 TR-Therapy device is equivalent to HCR 900 device, HCR 901, TECAR 901, MD-308.

Clinical Evaluation Report  
BTL-6000 TR-Therapy

We found out 5 articles published for the equivalent devices which confirm the intended use of BTL-6000 TR-Therapy. The articles were found and chosen within the literature search (*Comparative literature study for BTL-6000 TR-Therapy*) The approved intended use and intended population within these articles is summarized in the *Comparative literature study for BTL-6000 TR-Therapy* chapter 3.5, in the Table 3.

For all of these studies we found out whole articles, so the obtained clinical data are sufficient to prove the clinical performance and safety on. The abstract of articles is found out in the *Comparative literature study for BTL-6000 TR-Therapy* chapter 3.4.1.

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## Clinical Evaluation Report BTL-6000 TR-Therapy

### 5 Clinical safety

The clinical safety data was based on the search of the databases for adverse events reporting - MHRA and MAUDE and within the risk management process performed during the development for the BTL-6000 TR-Therapy device. The clinical safety data will be checked within the post market surveillance data after the product launch.

#### 5.1 Equivalent device

The data for the equivalent devices were collected within the literature search for the purpose of providing the evidence of their effective and safe clinical use – contact radiofrequency therapy.

Two (#18, #22) of the five chosen articles contain also the information to prove the clinical safety of the device. The first study proves the safety of the device based on the blood pressure measuring. The second study recorded the adverse events within whole study. Both articles confirm the safe use of the device. In first article (#18) any valuable deviations from the normal blood pressure weren't observed. In the second article (#22) any adverse events weren't reported within whole study.

We also checked within the literature search the databases for adverse events reporting - MHRA and MAUDE. Any reported adverse events weren't found.

#### 5.2 Post Market Surveillance Data

Activities of post market surveillance are set forth in *Post Market Surveillance Plan for the BTL-6000 TR-Therapy*. The results from the PMS will be described in the *PMS Report for BTL-6000 TR-Therapy*. It is compiled once a year. The report is designed usually at the beginning of the year in the January, it include post market surveillance data which are collecting during the year before. If some significant changes of the clinical data are discovered, the *Clinical Evaluation Report* is duly actualized. Because the BTL-6000 TR-Therapy the new technology in BTL, we declare the PMCF will be performed during next year for this device. The PMDF plan is described within the PMS activities (*Post market Surveillance Plan for BTL-6000 TR-Therapy*).

#### 5.3 Risk Analysis Review

The risk analysis identified risks associated with the use of the device. A multi-functional team considered the use of the device and identified the characteristics related to the safety of the medical device. The identification of hazards and the sequence of event(s) leading to hazardous situation(s) were conducted from both the product and user perspective (including foreseeable misuse). Both normal and fault conditions for the device were considered. Mitigation measures (controls) have been put in place to minimize the risk to the patient and user.

Risk management was created in November 2014. The risk were described and analyzed with FMEA method. The risk management include also the evaluation of the clinical risks.

##### 5.3.1 Residual risks

The device was designed with respect to the maximum safety of the operator as well as the patient. The BTL-6000 TR-Therapy device complies with all applicable safety standards.

The overall residual risk of the device is considered in accordance with the acceptability criteria as broadly acceptable. The evaluation of all risk management activities can declare that the health benefits outweigh the risks related to the device under the condition that the correct intended use and indications are followed. The risk analysis performed by the FMEA method confirmed that the device comply with all requirements for safety and essential performance.

Based on the assessment of all risk management activities it can be stated that medical benefits under the conditions of following the established intended use significantly outweigh possible risks related to the use of the device.

Risk analysis did not reveal any serious unacceptable risk concerning the product basic safety and essential performance. All risk management activities are in compliance with standard EN ISO 14971:2012.

The evaluation of the residual risk is described in the document *Risk management report* for BTL-049 (BTL-6000 Tr-Therapy).

231

## Clinical Evaluation Report BTL-6000 TR-Therapy

### 5.3.2 Clinical risks

The residual risks resulting from the risk management documentation concerning the clinical data are described in the user's manual. The clinical risks were assessed by the FMEA method within the risk management process. Following clinical risks were included for the evaluation:

IO1: Insufficient specification of intended use: The user's manual contains the information about the intended use as is approved via this clinical evaluation report.

IH04: Side effects and hazards are not specified: Among the clinical risks we have to state the side effects of tecar therapy. Within the literature search we found out in the competitor's documentation the possible side effects. These side effects are described in the user's manual for BTL-6000 TR-Therapy.

UE17: Not respecting the patient feedback during therapy. This risk was evaluated via the risk management process. The solution is the properly described handling with the device in the user's manual. This risk is described as a safety precaution in the user's manual.

UE29: Use of scratched capacitive electrode which could lead to local discomfort, pain. This risk was evaluated and solved via the risk management process with the internal testing and the safety precautions written in the user's manual.

UE30: Application of the therapy on the anesthetized patient or the patient with elevated pain threshold can cause local burns. This risk is described as a safety precaution in the user's manual. From the point of view of residual risk it is necessary to consider that the operator does not read the Manual or intentionally behaves in conflict with the User's Manual. In this case it is not possible to effectively eliminate the risk by safety precautions in the user's manual.

BI01, BI02: The biological hazards are solved via the risk management process. The materials which are in direct contact with the patients are in compliance with the standard ISO 10993-1. The biological tests were performed in the accredited laboratory.

TE01, TE09, TE10: Thermal risks that could cause the burn of patients: High temperature of applied parts during therapy, High temperature of surfaces during normal operation, Burn of the operator or patient. This is the common shown complaint for the thermal energy therapy devices. This risk was evaluated and solved via the risk management process with the internal testing. The internal performance tests prove the safety for this risk. Applied parts are not heated during the therapy.

The instructions in the user's manual are written in accordance with the clinical residual risks.

### 5.4 Conclusion – clinical safety

The clinical safety was supported by the risk management process for the BTL-6000 TR-Therapy and by the long term experiences with equivalent devices manufactured by Indiba supported also with the safety outcomes described in the chosen articles. We didn't find any adverse events or recalls for the equivalent device. The risk analysis performed within the development of the BTL-6000 TR-Therapy by the FMEA method confirmed that the device complies with all requirements for safety and essential performance. Based on all above written information we can state that the device is safe for use.

232

## Clinical Evaluation Report BTL-6000 TR-Therapy

### 6 Clinical performance

The clinical performance data was based on the long term experiences with equivalent devices manufactured by Indiba and within the documents regarding the manufacturing process for BTL-6000 TR-Therapy. These documents includes: the list of used components, applied parts of the device, materials used, internal and external testing, final inspection, development documentation and other steps of manufacturing process.

#### 6.1 Equivalent device

We chose 5 articles to prove the clinical efficiency for the intended use declared for BTL-6000 TR-Therapy. The articles were evaluated within the literature search. Only the articles which contain sufficient information to be able to undertake a rational and objective assessment were chosen. First article (#14) confirms the intended use for the Treatment of musculoskeletal acute and chronic injuries, acute and chronic diseases of musculoskeletal system. They studied patients most confirmed an improvement in pain symptoms. The second article (#15) deal with the Treatment for patients undergoing surgery of musculoskeletal system, post-surgical treatment fractures of the femur. The positive effect for pain and oedema in patients with femur fracture treated surgically was approved. The third article (#17) confirms the intended use such as Treatment of acute and chronic pathologies in sports, treatment of sports pathologies of bones, joints and muscles, Articular pathologies, muscular pathologies, tendinous pathologies. The fourth article (#18) shows the favourable results for the Treatment of degenerative neck pain. The fifth article speaks about the pain relieve in patients with isolated long head biceps (LHB) tendinopathy. The article shows mainly the effect of tecar in short-term follow up. All studies were performed for the adult patients.

Based on above written information we can declare the intended use for the BTL-6000 TR-Therapy as: treatment of musculoskeletal system acute, chronic and degenerative disorders and injuries such as: sports pathologies of bones, joints and muscles; tendinous, muscular and articular pathologies; neck pain; pain in patients with isolated long head biceps tendinopathy; post operative pain and oedema treatment of femur fractures.

#### 6.2 List of components

The list of all components used is described in the chapter 8 in the Technical file BTL-049. The material which is in the direct contact with the patient is described in the chapter 16 in the Technical file BTL-049. This chapter include also the biocompatibility test reports which were performed for the materials that are in the direct contact with the patients.

#### 6.3 List of internal and external tests

The list of all internal and external tests is described in the chapter 14 in the Technical file BTL-049.

The list of internal and external tests includes the internal performance tests, test issued by the accredited laboratories according to the applied safety standards.

The device was designed with respect to the maximum safety of the operator as well as the patient BTL-6000 TR-Therapy comply with all applicable safety standards. The safety of the device is verified by the testing performed according to the international electrical safety standard IEC 60601-1.

The materials in direct contact with the patients meet the requirements for the biocompatibility according to the ISO 10993-1, ISO 10993-5, ISO 10993-10. The sensitization, cytotoxicity and irritation test were performed for these materials. The test reports are stored in the chapter 16 in the Technical File BTL-049.

#### 6.4 Manufacturing process

Manufacturing Flow Chart provides a general understanding of the manufacturing process. The flowchart is stored in the chapter 11 in the Technical File BTL-049.

The company BTL Company is certified in accordance with EN ISO 13485:2012 with QMS certificate No. 75565-2010-AQ-CZS-SINCERT issued by DNV.

The company also meets the requirements for the Canadian quality management system. The company is certified with the QMS certificate No. CERT-0059956 issued by SAI GLOBAL.

The device has stored the documentation regarding the development of device in the Design History File for BTL-049.

133

## Clinical Evaluation Report

### BTL-6000 TR-Therapy

The final inspections for the devices are performed according to the final inspection templates – see chapter 9 in the Technical File BTL-049.

The risk analysis is described in the chapter 4 in the Technical File BTL-049.

Compliance with the essential requirements of MDD 93/42/EEC is described in the chapter 3 in the Technical File BTL-049.

#### 6.5 Conclusion – clinical performance

The clinical performance data was based on the long term experiences with equivalent devices manufactured by Indiba and within the documents regarding the manufacturing process for BTL-6000 TR-Therapy. The articles chosen within the literature search support us with the sufficient information to confirm the intended use of the BTL-6000 TR-Therapy as follows: treatment of musculoskeletal system acute, chronic and degenerative disorders and injuries such as: sports pathologies of bones, joints and muscles; tendinous, muscular and articular pathologies; neck pain; pain in patients with isolated long head biceps tendinopathy; post operative pain and oedema treatment of femur fractures. The device is intended for the adult patients.

The BTL Company has set the manufacturing process which meets the requirements of the Quality management standards in the Europe and in the Canada. We declare within the development documentation that the evaluated device is safe, suitable for intended use and complies with requirements defined in Annex 1 MDD 93/42/EEC (including article 6a).

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139

## Clinical Evaluation Report BTL-6000 TR-Therapy

### 7 User's manual

The device documentation include User's manual is in compliance with this clinical evaluation report.

#### 7.1 User's manual - Device intended use and contraindication

The users' manual contain following information regarding the clinical data:

The device is intended for treatment of musculoskeletal system acute, chronic and degenerative disorders and injuries such as: sports pathologies of bones, joints and muscles; tendinous, muscular and articular pathologies; neck pain; pain in patients with isolated long head biceps tendinopathy; post operative pain and oedema treatment of femur fractures.

The device is intended for the adult patients.

Contraindications for contact radiofrequency current: Patients with a pacemaker or other electronic implants must be approved beforehand by a specialist in a particular medical field. Patients must not touch the electrical conductors or grounded object, formation of undesirable pathways could occur for the transmission of high-frequency currents. People with low blood pressure or patients with treatment of the cervical vertebrae, spine or joints have to take a longer lasting therapeutic procedure, the therapy can lead to further reduction of blood pressure. Use of the device is contraindicated for pregnant women. When using resistance mode, it is necessary to remove the dentures if treatment is carried out in the relevant area. Patients or the users should remove rings, necklaces, piercings and other metal objects. Do not apply on damaged skin, for example in the case of open wounds or fresh burns.

The intended use is based on the particular informations found in the chosen articles for the equivalent devices. The contraindications and side effects are based on the informations which the competitors declare within the labelling for the equivalent tecar therapy devices. We collected information from all 4 equivalent devices and based on them we set the contraindications and side effect for BTL-6000 TR-Therapy device.

#### 7.2 User's manual - Warnings

The warning and instruction in the user's manual are in compliance with the residual risks and clinical risks obtained within the risk management documentation and clinical evaluation. Other warnings in the manual arise from foreign countries legislation requirements or from the customer requirements. The user's manual contains the clinical data evaluated within this clinical evaluation report.

#### 7.3 Conclusion – User's manual

The user's manual contains the same intended use of the device and contraindications as is supported by this clinical evaluation report. The user's manual contains enough information regarding the clinical performance and safety of the device.

## Clinical Evaluation Report BTL-6000 TR-Therapy

### 8 Conclusion

#### 8.1 Intended Use Benefit evaluation

With the increasing awareness of the dangers of the side-effects of pharmaceutical therapies, contact radiofrequency therapy is starting to get attention as one of the forms of alternative treatment. Chiropractors and physical therapists use contact radiofrequency therapy to improve the quality of life through pain management and biostimulation. Contact radiofrequency therapy can also reduce the tightening of muscles which are often responsible for significant pain in the back. Reducing muscle spasms and pain in your back can prevent exacerbation of the back problems, which may lead to disability or surgical intervention.

General benefits resulting from the above mentioned truths of BTL-6000 TR-Therapy are mainly the non invasive treatments of a many different dysfunctions. The main use is focused on the movement or nervous dysfunction. The patient's recovering is faster when is supported by this treatment.

#### 8.2 Conclusion - Clinical risk / benefit

The overall residual risk of the device is considered in accordance with the acceptability criteria as broadly acceptable. The evaluation of all risk management activities can declare that the health benefits outweigh the risks related to the device under the condition that the correct intended use and indications are followed. The risk analysis performed by the FMEA method confirmed that the device comply with all requirements for safety and essential performance. Based on the assessment of all risk management activities it can be stated that medical benefits under the conditions of following the established intended use significantly outweigh possible risks related to the use of the device.

#### 8.3 Overall evaluation

The clinical data obtained and experiences with the used technology for the equivalent devices were assessed as a sufficient for the confirmation of the performance and safety of the device.

The clinical safety was supported by the risk management process for the BTL-6000 TR-Therapy and by the long term experiences with equivalent devices manufactured by Indiba supported also with the safety outcomes described in the chosen articles. We didn't find any adverse events or recalls for the equivalent device. The risk analysis performed within the development of the BTL-6000 TR-Therapy by the FMEA method confirmed that the device comply with all requirements for safety and essential performance. Based on above written information we can state that the device is safety for use.

The clinical performance data was based on the long term experiences with equivalent devices manufactured by Indiba and within the documents regarding the manufacturing process for BTL-6000 TR-Therapy. The articles chosen within the literature search support us with the sufficient information to confirm the intended use of the BTL-6000 TR-Therapy as follows: treatment of musculoskeletal system acute, chronic and degenerative disorders and injuries such as: sports pathologies of bones, joints and muscles; tendinous, muscular and articular pathologies; neck pain; pain in patients with isolated long head biceps tendinopathy; post operative pain and oedema treatment of femur fractures. The device is intended for the adult patients.

The user's manual contains the same intended use of the device and contraindications as is supported by this clinical evaluation report. The user's manual contains enough information regarding the clinical performance and safety of the device.

The intended use is based on the particular informations found in the chosen articles for the equivalent devices. The contraindications and side effects are based on the informations which the copatitors declare in the user's manuals for the equivalent tecar therapy devices.

Indeed, this fact allows us to conclude that the BTL-6000 TR-Therapy devices intended for performing physiotherapy by contact radiofrequency therapy are suitable for their state-of-the-art intended use, to be marketed as BTL-6000 TR-Therapy, provided that the requirements for their specified characteristics and performance, under the intended conditions and for its intended purpose of use, are fulfilled. Summarized, BTL-6000 TR-Therapy does not pose any undue safety concerns either to the operator, physician or patient and the risks associated with the use of this device are acceptable when weighed against the benefits to the patient.

236

Clinical Evaluation Report  
BTL-6000 TR-Therapy

**9 Documentation – cross reference**

**9.1 Unpublished literature**

- Complaints and Adverse Events electronic databases:
- MHRA (<http://www.mhra.gov.uk/>)
- MAUDE (<http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfmaude/search.cfm>)
- The websites of the competitive producers:
- Indiba (<http://www.indiba.com/>)
- Human Tecar (<http://www.humantecar.eu/?lang=en>)

**10 Statement**

We declare within this unpublished/published documentation that product BTL-6000 TR-Therapy is safe and suitable for intended use. In accordance with the definition of clinical data set forth in Medical Devices Directive 93/42/EEC, as amended by 2007/47/EC, we have chosen the way of collecting clinical information about safety and performance following from a device's intended use: Through published and/or unpublished written documents concerning clinical experience with equivalent devices. The analogy/equivalence of the BTL-6000 TR-Therapy device to currently used equivalent devices was proved via the clinical, biological, technical equivalence. The device complies with requirements defined in Annex 1 MDD 93/42/EEC article 1, 3 and 6a.

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© Все права защищены. Никакая часть данного руководства не может быть воспроизведена, сохранена в научно-исследовательском центре или передана любыми средствами, включая электронные, механические, фотографические или другие документы, без предварительного согласия с BTL Industries Limited.